

All-Ceramic Chairside Preparation Guide for IPS e.max®

For IPS e.max® indications

		IPS e.max Ceram/ZirPress	IPS e.max ZirCAD	IPS e.max CAD/Press
Anterior	Crown		✓	V
	Veneer	✓	✓	V
	Bridge		V	V
Posterior	Crown		V	V
	Bridge		✓	


Anterior Chairside Preparation Guide


Full-Coverage Restorations

ANTERIOR CROWN PREPARATION

3-Unit Bridge Restorations

3-UNIT BRIDGE PREPARATION





Veneers

A medium grit, round-ended, diamond bur is used to remove a uniform thickness of facial enamel by joining the depth-cut grooves.

The diamond bur is angled to bevel back the incisal edge.



Thin Veneers

IPS e.max can be pressed to as thin as 0.3 mm for veneers. If sufficient space is present, IPS e.max can be placed over the existing teeth without the removal of any tooth structure. Depending on the case requirements, however, some teeth may need to be prepared to accomodate for the thickness of the ceramic and to ensure for proper contour and emergence profile.





All-Ceramic Chairside Preparation Guide for IPS e.max®

Posterior Chairside Preparation Guide

Inlays/Onlays


INLAY PREPARATION

1.5 mm to 2.0 mm isthmus width
Shoulder margin

Rounded internal
line angles

1.0 mm- to 1.5 mmwide gingival floor

1.5 mm-depth
at isthmus



3-Unit Bridge Restorations

3-UNIT BRIDGE PREPARATION



1.5 mm axial reduction

Rounded internal line angles Shoulder margin 1.5 mm axial reduction 1.0 mm reduction at the gingival margin Shoulder margin

Full-Coverage Restorations

POSTERIOR CROWN PREPARATION

Rounded internal line angles

1.5 mm to 2.0 mm occlusal reduction

1.5 mm axial reduction

1.0 mm reduction at the gingival margin margin A flat-ended, tapered diamond is utilized to establish a shoulder margin

When layered or pressed ceramic margins are preferred in conjuction with a zirconia framework, enhanced gingival esthetics can be achieved with a definitive 90 degree shoulder preparation.

www.ivoclarvivadent.com

Call us toll free at 1-800-533-6825 in the U.S., 1-800-263-8182 in Canada.

©2009 Ivoclar Vivadent, Inc. IPS e.max is a registered trademark of Ivoclar Vivadent.

#598204 All illustrations ©2009 Montage Media Corporation

CONVENTIONAL CEMENTATION PREPARATION

Occlusal reduction of at least 2.0 mm in contact area

Shoulder preparation of at least | 1.0 mm proximal reduction | Coronal length at least 4.0 mm

