

# Six Levels Of Protection



**Junior:** 1 layer of EVA material (3mm) with added incisal and occlusal protection. Designed specifically for children with mixed dentition.



**Light:** 2 layers of laminated EVA (3mm) material. Designed for wrestling, volleyball, mountain biking, and motocross. *(Only available in clear, bright yellow, bright red, light blue or green.)*



**Light Pro:** 3 layers of laminated EVA (3.5mm) material (one hardened layer). Specially designed to provide maximum comfort and protection for all sports.



**Medium:** 2 layers of laminated EVA (5mm) material. Specially fabricated for soccer, rugby, basketball, softball, rollerblading, and skateboarding.



**Heavy:** 2 layers of laminated EVA (5mm) material with 3 unique power dispersion bands. Specially designed for baseball, football, racquetball, martial arts and boxing.



**Heavy Pro:** 3 layers of laminated EVA (5mm) material (one hardened layer). Custom fabricated for ice, field and street hockey along with kick-boxing & other heavy contact sports where blows from pointed objects are expected.

•*Helmet straps available in red and black*